

LABORHITS

SGE's GC Inlet Liner Range



Easy to Choose

Color coded by geometry to simplify your selection.

Confidence in Your Analysis

Certified deactivation gives you confidence in your analysis.

Easy to Use

Contain o-rings so you're ready to go.

Spot the Difference in the SGE Lineup!

Choosing the right inlet liner and injection parameter can increase peak areas and reduce detection limits by up to 300%.

Customer research found that a significant number of GC users don't understand the importance of inlet liner selection, or how it contributes to their analysis. The SGE inlet liner range aims to make it simple for all gas chromatographers to select the right liner.

To select the right inlet liner for your instrument, injection and sample use the SGE liner selection tool: www.sge.com/linertool

